

Selectin Biosciences Inc.

β (1-4)-Galactosidase

Recombinant from *Streptococcus pneumonia*

β (1-4) Galactosidase (β -D-galactoside galactohydrolase, EC 3.2.1.23) releases only β (1-4)- linked, non-reducing terminal galactose from complex carbohydrates and glycoproteins. β (1-4) galactose is by far the most common linkage found in N-linked oligosaccharides. The enzyme is as active on tetraantennary oligosaccharides as on disaccharides containing β (1-4)-linked galactose. Fucose linked to the penultimate N-acetylglucosamine will block cleavage of the galactose. Up to 100 μ g of asialofetuin can be completely β (1-4)-degalactosylated in less than 1 hour using 3 mU of enzyme.

β (1-4) Galactosidase is isolated from a clone of *Streptococcus pneumonia* expressed in *Escherichia coli*.

β (1-4) Galactosidase is useful for:

- Structural analysis of oligosaccharides
- Distinguishing different galactose linkages
- Removing heterogeneity from glycoproteins

Product Code: GE 12

Specifications

Activity: ≥ 6 U/mg, ~ 3 U/mL

Storage: Store at 4°C.

Formulation: The enzyme is provided as a sterile solution in 20 mM Tris pH 7.5, 25 mM NaCl.

Stability: Stable at least 12 months when stored properly. Several days exposure to ambient temperatures will not reduce activity.

Product Description

Molecular Weight: 240,000 Daltons

Purity: Each lot of β (1-4) Galactosidase is tested for contaminating activities by incubating the enzyme for 24 hours at 37°C with 2 μ L of enzyme. SDS-PAGE analysis of the treated BSA shows no evidence of degradation.

The production host strain has been extensively tested and does not produce any detectable glycosidases.

Specificity: Non-reducing terminal β (1-4) galactose. Number of antennae does not affect cleavage rate. Fucose linked to the penultimate N-acetylglucosamine will block cleavage of the galactose.

Assay

One unit of β (1-4) Galactosidase is defined as the amount of enzyme required to produce 1 μ mole of p-nitrophenol (pNP) in 1 minute at 37°C, pH 5.0 from p-nitrophenyl- β -D-galactopyranoside.

Reagents

- 5X Reaction buffer 6.0 – 250 mM NaHPO₄, pH 6.0

Selectin Biosciences Inc.

Suggestions for Use

Procedure For Degalactosylation

1. Add up to 100 μ g of asialoglycoprotein or 1 nmol of oligosaccharide to tube.
2. Add water to 14 μ L
3. Add 4 μ L 5X Reaction Buffer
4. Add 2 μ L β (1-4) Galactosidase
5. Incubate at 37°C for 1 hour.

For glycoproteins, cleavage may be monitored by SDS-PAGE if the size differential between native and de-galactosylated protein is sufficient for detection.

Note: The optimum for cleavage of oligosaccharides is ~ pH 6.0.

References

1. Glasgow, L.R., J.C. Paulson and R.L. Hill. Systematic purification of five glycosidases from *Streptococcus pneumoniae*. **J. Biol Chem** **252**: 8615-8623(1977).
2. Kobata, A. Use of endo- and exoglycosidases for structural studies of glycoconjugates. **Anal Biochem** **100**: 1-14 (1979).
3. Prime, S. J. Dearnley, A.M. Venton, R.B. Parekh and C.J. Edge. Oligosaccharide sequencing based on exo- and endoglycosidase digestion and liquid chromatographic analysis of the products. **J Chromatogr A** **720**: 263-274 (1996)

4. Dwek, R.A. , C.J. Edge, D.J. Harvey, M.R. Wormald and R.B. Parekh. Analysis of glycoprotein-associated oligosaccharides. **Ann Rev Biochem** **62**: 65-100.

Selectin Biosciences Inc warrants that the above product conforms to the specifications described herein. Should the product fail for reasons other than through misuse Selectin Biosciences Inc. will, at its option, replace free of charge or refund the purchase price. This warranty is exclusive and Selectin Biosciences Inc. makes no other warrants, expressed or implied, including any implied conditions or warranties of merchantability or fitness for any particular purpose. Selectin Biosciences Inc. shall not be liable for any incidental, consequential or contingent damages.

This product is intended for in vitro research only.

REVISION 5/22/12